

Instructions for Use

Soleo*line*


Soleo **SonoStim**, Soleo **Galva**, Soleo **VacoS**



Explanation of Symbols

Soleo SonoStim / Soleo Galva
Front panel

Fig. 1



Device and operating elements


- 1 Intensity controller channel I
- 2 Intensity controller channel II
- 3 Alligator clips clipping options
- 4 Tray mount for ultrasound head
- 5 Display
- 6 Touch pen in holder
- 7 Slot for SD card

Explanation of Symbols

Soleo SonoStim / Soleo Galva

Display/ Navigation bars


Fig. 2



Displays

- 8 Status bar
- 9 Buttons on the screen
- 10 Title bar
- 11 Navigation bars

Fig. 3



Navigation bar

Description of functions


- | | |
|----------------------|---------------------------------|
| (A) Back | Moves back one step |
| (B) Memory | Switches to the memory area |
| (C) Indications | Switches to the indication menu |
| (D) Programmes | Switches to the programme list |
| (E) Favourites | Switches to the favourites area |
| (F) Scroll backwards | Moves back one page |
| (G) Scroll forwards | Moves forwards one page |

Explanation of Symbols

Soleo SonoStim / Soleo Galva

Rear Panel/Switches and Connector Sockets

Fig. 4



Switches and
connector sockets

- 12 On/off switch
- 13 Holder for mains fuse
- 14 Connection for mains cable
- 16 Socket for electrode cable channel II
- 17 Connector for ultrasound head 0.8 / 2.4 MHz
- 18 Socket for electrode cable channel I


Note: The socket (15) has no use for the Soleoline.*

Explanation of Symbols

VacoS

Front panel/sockets

Fig. 5



Ports

Group 1

- 19 Socket for Vaco electrode tube black, cathode
- 20 Socket for Vaco electrode tube red, anode

Group 2


- 21 Socket for Vaco electrode tube black, cathode
- 22 Socket for Vaco electrode tube red, anode

Explanation of Symbols

VacoS

Rear Panel/Switches and Connector Sockets

Fig. 7




Switches and connector sockets

- 23 On/off switch
- 24 Connection for short mains connection cable
- 25 Connection for mains cable
- 26 Bracket for mains backup
- 27 Socket for electrode cable channel II
- 28 Socket for connection cable with SonoStim/Galva channel II
- 29 Socket for electrode cable channel I
- 30 Socket for connection cable with SonoStim/Galva channel I
- 31 Soleo SonoStim / Galva identification plate
- 32 Soleo VacoS identification plate

Explanation of Symbols

Soleo SonoStim / Soleo Galva / VacoS
Side view

Fig. 8



Connections


33 Water separator connector valve

Explanation of Symbols


Soleo SonoStim / Soleo Galva
Ultrasound heads

Fig. 9

Ultrasound head, large



Ultrasound head, small



Explanation of Symbols



In the instructions for use, this symbol indicates "Danger".

Caution!

In the instructions for use, this symbol represents "Caution" with regard to damage of the device.



Device Type BF



Follow the instructions for use



Instructions for Use



Serial number



Item number




Manufacturer



Date of manufacture



This symbol indicates danger areas on the machine.



Interval operation

Contents

Explanation of Symbols

Soleo Sono*Stim*/Soleo *Galva*

Front panel

Display/ Navigation bars

Rear panel / Switches and connector sockets

VacoS

Front panel/sockets

Rear panel / Switches and connector sockets

Soleo Sono*Stim*/Soleo *Galva* and VacoS

Side view

Ultrasound heads

Explanation of Symbols

| | | Page |
|----------|---|------|
| 1 | Indications/Contraindications | |
| | 1.1 Electrotherapy indications | 1 |
| | 1.2 Electrotherapy contraindications | 3 |
| | 1.3 Ultrasound therapy indications | 4 |
| | 1.4 Ultrasound therapy contraindications | 5 |
| 2 | Side Effects | 6 |
| 3 | Application Information | |
| | 3.1 General | 7 |
| | 3.2 Electrotherapy | 8 |
| | 3.3 Ultrasound | 9 |
| 4 | Warnings | |
| | 4.1 General | 10 |
| | 4.2 Electrotherapy | 11 |
| 5 | Soleoline – in brief | 12 |
| 6 | Intended Use | 13 |
| 7 | 7 Device Set up | |
| | 7.1 Installation of Soleo Sono <i>Stim</i> /Soleo <i>Galva</i> cables | 14 |
| | 7.2 Installation of the VacoS cable | 15 |

Contents

| | | |
|-----------|--|--|
| 8 | Basic Settings 8.1 General 8.2 Electrotherapy 8.3 Ultrasound Therapy 8.4 Maintenance | 17 21 22 23 |
| 9 | Soleo SonoStim - Instructions for operation 9.1 Electrotherapy 9.2 Ultrasound Therapy 9.3 Water Bath Treatment 9.4 Simultaneous Therapy | 24 29 33 34 |
| 10 | Soleo Galva - Instructions for operation 10.1 Electrotherapy | 37 |
| 11 | Soleo SonoStim/Soleo Galva - Instructions for operation 11.1 Electrotherapy with VacoS 11.2 Water separator VacoS 11.3 SD card 11.4 Menu indications 11.5 Storing Programmes 11.6 Retrieving programmes/Editing the list 11.7 Programme sequence | 38 41 42 43 45 47 49 |
| 12 | Paralysis Diagnosis and Therapy 12.1 Overview 12.2 Medium frequency test by Lange 12.3 Neofaradic test 12.4 Chronaxie/accomodation quotient 12.5 Fisch Gold Test | 51 52 53 54 57 |
| 13 | Technical Information 13.1 Soleo SonoStim/Soleo Galva - General 13.2 Stimulation 13.3 Ultrasound 13.4 VacoS | 59 60 62 63 |
| 14 | Cleaning / Disinfection | 64 |
| 15 | Electrodes - Information on use of the various types of electrodes | 66 |
| 16 | CE mark / manufacturer | 67 |

Contents

| | | |
|-----------|---|----------------|
| 17 | Scope of Delivery, Accessories 17.1 Soleo Sono <i>Stim</i> 17.2 Soleo <i>Galva</i> 17.3 VacoS | 68 69 70 |
| 18 | Device Combinations | 71 |
| 19 | Safety and Maintenance | 72 |
| 20 | Function Test | 73 |
| 21 | Legal Notices | 74 |
| 22 | Error Message, Troubleshooting, Disposal | 75 |
| 23 | Manufacturer's EMC Declaration | 79 |

Valid for devices Soleo Sono*Stim*, Soleo *Galva* and VacoS.

These instructions for use are an integral part of the device.
It must be stored with the device and kept accessible at all times for anyone to operate this device.

The instructions for use are valid from August 2016.

1.1 Electrotherapy indications

Circulation promotion treatment

- Functional circulatory disorders
- Angiopathy in diabetes mellitus
- Angioneuropathy (M. Raynaud)
- Acrocyanosis
- Arterial occlusive disease (stage I/IIa)
- Venous vascular disease
- Chronic venous insufficiency
- Sympathetic reflex dystrophy
- Sudeck syndrome (stage II)
- Disruption of trophism
- Lymphedema
- Fascia treatment (as pretreatment)
- Chronically tense muscles
- Disease of the peripheral joints (chronic)
- Radicular syndromes with spinal disorders (subacute and chronic)
- Ligament ruptures
- Tenosynovitis
- Bursitis

Analgesic therapy

- Sprains, contusions, muscle sprains,
- Tendons, joints and bones
- Large myalgia
- Polyarthralgia
- Joint diseases, inflammatory
- Chronic polyarthritis (rheumatoid arthritis), treatment in inflammatory poor intervals
- Diffuse pain in periarthropathy, arthrosis,
- M. Bechterew
- Tendinitis, tenosynovitis
- Local and pseudoradicular spinal disorders
- Radical lesions
- Disc surgery
- Chronic inflammatory spinal diseases
- Chronic pain conditions, pain disorder
- Pelvipathy (so-called adnexitis)
- Pelvipathy (so-called prostatitis)
- Neuralgia
- Neuritis
- Polyneuropathy
- Radicular syndrome
- Nerve compression syndrome
- Allodynia
- Causalgia

1.1 Electrotherapy indications

| | |
|-------------------------------------|--|
| Trophism promotion treatment | <ul style="list-style-type: none"> • Venous ulcers • Neurotrophic ulcers • Decubitus ulcers • Delayed wound healing • Osteoporosis • Delayed bone healing • Ligament ruptures |
| Muscle toning therapy | <ul style="list-style-type: none"> • Muscle weakness • Voluntary innervation weakness • Missing muscle feeling |
| Muscle detoning therapy | <ul style="list-style-type: none"> • Overloading the muscles • Painful muscular tension • Myofascial pain syndromes • Tendon and muscle rupture • Contractures |
| Other indications | <ul style="list-style-type: none"> • Periarthropathies and insertion tendinopathies (acute, subacute) • Periarthropathies and insertion tendinopathies (chronic) • Disease in the peripheral joints (acute, subacute) • Post-traumatic state images (hematoma, edema) • Excessive sweating of hands and feet (hyperhidrosis palmarum and plantarum) • Spinal disorders acquired: post-traumatic, post inflammatory, degenerative • Posture spinal disorders in childhood and adolescence • Spinal fractures, spinal fusion • Joint surgery, arthroscopy, arthrotomy • Congenital malformations or incorrect settings of the musculoskeletal system to be treated in childhood • Fecal incontinence • Urinary incontinence • Disruption of gut motility • Chronic constipation • Secondary lymphedema and lymphoedema • Pelvic pain • Polyneuropathy • Anterior horn of the spinal cord disease • Peripheral paralysis, brachial plexus • Non-local, generalised spinal disorders • Spasticity • Spastic paresis • Flaccid paralysis |

1.2 Contraindications for electrotherapy

General Contraindications

- Vague pain symptoms
- Acute inflammation (local, systemic)
- Suspected diseases of the cardiovascular system
- (Suspected) epilepsy
- Arterial occlusive disease from stage IIb according to Fontaine
- Purulent processes
- Fever
- General infections
- Tumours of a malignant and benign nature
- Thrombophlebitis
- DVT at risk of embolism
- Pacemakers and other implanted electronic devices
- Transcardial current flow
- Metal implants in the current field when using galvanic or unipolar currents with pulse widths of more than 1 ms
- Psychosis
- Reduced skin sensitivity
- Application on or in the vicinity of open wounds or fractures

The following should also be considered when planning TENS therapy

- Causally repairable pain
- Predominantly psychogenic pain (TENS therapy is ineffective)
- Central pain syndromes, for example thalamic pain syndrome

In addition the following contraindications should be observed

Muscle treatment with swollen group pulses (threshold currents) in case of:

- Unconsciousness, consciousness
- Reflex inhibition (e.g. fractures)
- Muscle inflammation
- Mimetic facial muscles
- If the electric stimulation therapy triggers persistent pain
- Infants, toddlers
- When spasticity increases
- Pregnancy

Particular caution is required:

- With bleeding after acute trauma or fracture
- After surgical procedures when muscle contraction may disrupt the healing process

Precautions for certain current forms

Diadynamic current, trabert current

- A careful electrode technique is necessary to avoid possible skin damage through a high galvanic share of diadynamic currents
- Cautious dosing in case of sensory disturbances

Electroacupuncture

- No needle acupuncture in bleeding disorders

1.3 Indications for ultrasound therapy

Indications from orthopaedics, surgery, traumatology, rheumatology

- Spinal pain syndromes, e.g. cervical spine syndrome
- Ankylosing spondylitis (only in inflammation-free intervals)
- Joint disorders
- Rheumatoid arthritis (if heat treatment is indicated)
- Arthrosis
- Periarthropathies
- Epicondylopathy
- Tendopathies, periostosis, calcaneal spurs
- Achillodynia
- Scars, contractures, Dupuytren's contracture
- Posttraumatic arthritis
- Fractures (particularly in cases of delayed callus formation)

Other indications

- Bronchial asthma
- Rhinopathy
- Persistent symptoms of the cervical spine after whiplash injury with recurrent blockages
- Headache
- Earache
- Post herpetic neuralgia
- Functional disorders of the stomach and duodenum
- Pelvic pain
- Functional disorders of the pelvis

1.4 Contraindications for ultrasound therapy

General Contraindications

- Vague pain symptoms
- Diseases where heat should not be used, for example, acute inflammatory diseases
- Diseases where mechanical impacts are contraindicated, e.g. deep vein thrombosis,
- Restricted blood flow
- Suspected diseases of the cardiovascular system
- Haemorrhagic diatheses
- Cervical vertebrae – do not use ultrasound higher than C 3
- Ultrasound treatment of parenchymatous or heat-sensitive organs (testes, eyes, gravid uterus, liver, kidneys etc.)
- Anaesthetised skin areas
- Temperature sensitivity disorders
- After treatment with ionising radiation
- Epiphyseal plates
- Tumours
- Electronic pacemakers
- Application on or in the vicinity of open wounds
- Impaired reflexes or pain sensitivity
- Application to the abdomen, pelvis, or lumbar area when pregnant or presumed pregnant

Metal implants and endoprotheses

Dynamic ultrasound radiation in low dosage is nowadays no longer a matter for concern.

**Side effects of
electrotherapy**

If used correctly, no side effects are known.

**Side effects of
ultrasound therapy**

If used correctly, no side effects are known.

3.1 General

Prior to the application of the device to the patient, the user should become acquainted with the operating instructions or the individual treatment methods and also with the indications / contraindications, warnings and directions for use. Additional sources of information on ultrasound therapy should also be taken into consideration.

This user manual should always be kept with the device to enable the persons charged with the operation of the device to access it at any time.

Caution! To avoid the risk of electric shock, the device must be disconnected from the mains supply before performing any maintenance or cleaning activities.

Caution! Before use, ensure that the device is powered via a properly earthed plug with a grounded outlet (electrical installation according to DIN VDE 0100 Part 710). The device must only be operated with the supplied power cord. The power cord must be protected against mechanical stress.

Caution! Magnetic and electric fields can influence the function of the device. Therefore, the device must not be operated in the vicinity of devices that generate strong electromagnetic fields (X-ray or diathermy devices, MRI scanners). Please keep a safe distance of several metres.

Caution! The device is not suitable for use in areas with an explosive, combustible or combustion-promoting atmosphere.

Caution! During use, the device must be placed in such a way that direct access to the device's central power supply is possible, so that it can be disconnected from the mains at any time.
Inspect the device before use. If damaged, it must not be used.

Caution! Only accessories provided by Zimmer MedizinSysteme GmbH are to be used and no other devices must be connected to the VacoS.

Caution! The device may cause malfunctions or may interfere with the operation of equipment in the vicinity. It may be necessary to take appropriate remedial action, such as the realignment or re-arrangement of the device or the shield.

Caution! No change to the device or the medical system may be carried out during the product life of the device.

The entire medical system is suitable for application in the patient environment.

Caution! To safely disconnect the unit from the mains supply, remove the mains plug from the socket housing or socket.

3.2 Electrotherapy

Electrode positioning and choice of treatment parameters should always be adapted to the treatment guidelines.

For currents in which a chemical burn is possible (e.g. galvanic current, diadynamic current, currents with galvanic share), the recommended maximum current density is 2 mA rms/cm² electrode surface.

If the current density exceeds 2 mA eff/cm², an increased attention from the user is required.

Note that with a smaller electrode surface a higher current density is achieved by using different electrodes.

The Tens treatment should not be used on the head or directly on the eye.

Caution!

In simultaneous operation, the ultrasonic head may be used as the dynamic electrode in connection with a stationary electrode, which is connected to the electrode cable. The ultrasonic and electrode channels form an application part. Please note that therapeutic currents flow in this mode via the ultrasound head.

Caution!

3.3 Ultrasound

Handle the ultrasound head carefully, as rough handling can alter its properties. Do not bring the ultrasound head into contact with sharp or pointed objects, as the aluminium head is susceptible to scratching.

The use of a couplant other than the special Sono plus ultrasound gel can damage the ultrasound head.

Disinfect the ultrasound head after use, using commercially available equipment disinfectants.

Inspect the ultrasound head for damage prior to each use. If damaged, it must not be used.

Caution!

In simultaneous operation, the ultrasonic head may be used as the dynamic electrode in connection with a stationary electrode, which is connected to the electrode cable. The ultrasonic and electrode channels form an application part. Please note that therapeutic currents flow in this mode via the ultrasound head.

Caution!

With regular use of the ultrasound head in the water bath, the user must take protection measures. All body parts of the user immersed in the water must be protected by appropriate measures. Suitable protective measures are so-called "aerated clothes".

This can be neoprene gloves or latex over wool gloves.

4.1 General



This product is intended for use by medical professionals only.



Very important: Never connect two patients to the device in one treatment session!



The patient must not be left unattended during therapy.



Treatment instructions regarding the treatment location, duration and intensity of treatment require medical knowledge and may only be issued by licensed physicians, therapists and members of the auxiliary medical professions. It is imperative that these instructions are followed.



Use in wet areas is not permitted and non-compliance may lead to considerable damage and endanger both the patient and the user.



Dispose of the packaging material correctly. Keep out of reach of children.



The use of this product with settings or for applications other than those specified in the user manual may lead to hazards resulting from the uncontrolled release of ultrasonic energy.



In the case of patients with implants or with an implanted electronic device, treatment may only be performed after clarification that no risk is involved.



The simultaneous connection of the patient to a high-frequency surgical device is not permitted. This can lead to burns.

4.2 Electrotherapy

Electrotherapy

Use of electrodes near the thorax may increase the risk of ventricular fibrillation.

When performing an iontophoresis, the drug used may possibly have an analgesic effect and pain sensitivity is then reduced.

Stimulation current can have a stimulating effect on insulin release. In diabetes patients this may therefore lead to hypoglycaemia.

The electrical stimulation or the materials used can cause skin irritation or hypersensitivity in susceptible patients. This can be reduced by the use of alternative electrode material or by changing the electrode system.

Stimulation should not be used:

- above the sinus nerve
- on the throat and mouth
- transthoracically
- transcerebrally
- over swollen, infected, inflamed areas
- over or in cancerous lesions

Vacuum

In patients who are prone to bruising, perform an evaluation of the lack of risks before beginning treatment.

Note:

For electrical therapy in combination with the vacuum unit, the above warnings apply,



With an intensity regulator, currents above 10 mA eff. or voltages in excess of 10 V at the output terminals are permitted.

| | |
|-----------------------------------|--|
| What is Soleoline? | <p>A ultra-modern innovative product line that provides 3 different device variants.</p> <p>Soleo SonoStim A state of the art, innovative combination device for electrotherapy and ultrasound therapy, with the option of the additional use of a vacuum unit.</p> <p>Soleo Galva A state of the art, innovative electrotherapy device, with the option of the additional use of a vacuum unit.</p> <p>Soleo Sono A cutting-edge, innovative ultrasound therapy device.</p> <p><i>Note:</i> <i>The operation of Soleo Sono is described in a separate manual.</i></p> |
| What does Soleoline offer? | <p>The delivery of monophasic, biphasic and medium frequency currents for nerve stimulation and muscle therapy in 1-channel and 2-channel operation and the delivery of therapeutic ultrasound.</p> |
| Why use Soleoline? | <p>A modern, clearly arranged colour display that provides all therapy-related parameters, as well as modern touch operation. Customisable program start-up settings and a clear, simple menu offer maximum convenience for the user.</p> <p>The combination of electrotherapy and ultrasound therapy in one system allows proven simultaneous therapy.</p> <p>The compact design allows space-saving work in the practice room and is ideally suited for use during home visits.</p> <p>The use of the vacuum unit ensures comfortable electrode application and additionally results in a pleasant massage effect for the patient.</p> |
| Innovations in Soleoline? | <p>SonoSwing, innovation in the field of ultrasound therapy:</p> <ul style="list-style-type: none"> - an ultrasound head with two frequencies - selectable penetration depth by the percentage adjustment of the frequency components. <p><i>Note:</i> <i>The application of the device is intended for the medical community (e.g. doctors, therapists and trained medical assistants).</i></p> |

The *Soleoline* product range includes 3 different device versions:

Soleo *Sono*

Ultrasound therapy device for treatment with therapeutic ultrasound.

Soleo *SonoStim*

Combination device for ultrasonic therapeutic treatment and electrotherapy, with the option of the additional use of a vacuum unit.

Soleo *Galva*

A device for electrotherapy, with the option of the additional use of a vacuum unit.

7 Device Set up

7.1 Installation of cables

Soleo SonoStim / Soleo Galva

7

Note: *The following description of the installation of the cables refers to the operation of Soleo SonoStim/Soleo Galva without VacoS.*

Note: *There are coloured arrows located on the cables as a guide for proper connection to the device.*

Electrotherapy When connecting the electrode cable, ensure that the green arrow is in the plug below.

- Connect the electrode cable for channel I to the appropriate socket (18).
- Connect the electrode cable for channel II to the appropriate socket (16).
- Insert the red alligator clips to the red connector of the electrode cable.
- Insert the black alligator clips to the black connector of the electrode cable.

Ultrasound therapy Connect the ultrasound head to the appropriate socket (17).

Connect the mains cable Connect the mains cable to the appropriate socket (14) and connect the cable to the mains.

Note: *The device may only be connected to power outlets with a protective contact.*

Switch the device on The device is switched on via the rocker switch (12).

Switch the device off The device is switched off via the toggle switch (12).
To completely disconnect the device (all-phase) from the network, remove the power cord.

Caution! All cables must be protected against jamming or other mechanical damage.

Note: *All of the buttons, menus and sub-menus can be activated directly on the screen using your finger or a touch pen.*

7.2 Installation of cables for VacoS

| | |
|--|---|
| Note: | <i>The control and operation of VacoS is via Soleo SonoStim/Soleo Galva. A stand-alone operation is not possible.</i> |
| Preparation | Position Soleo SonoStim/Soleo Galva so that the devices are aligned with each other. The sockets of the Vaco electrodes hoses must be on the front (fig. 4). |
| Connection of VacoS with SonoStim/Galva | <p>Make sure when mounting the connecting cable that when connecting the cable to Soleo SonoStim/Soleo Galva (16+18), the green arrow is in the plug below.</p> <p>Make sure when mounting the connecting cable that when connecting the cable to VacoS (28+30), the green arrow is in the plug below.</p> <p>Channel I Connect the connecting cable to the intended port on Soleo SonoStim / Soleo Galva (18) and on VacoS (30).</p> <p>Channel II Connect the connecting cable to the intended port on Soleo SonoStim / Soleo Galva (16) and on VacoS (28).</p> |
| Note: | <i>If Soleo SonoStim / Soleo Galva is operated in conjunction with VacoS, the connection is made with the electrode cable on VacoS. The assembly is described below.</i> |
| Assembly Electrode cable | <p>Make sure when mounting the electrode cable for VacoS that the green arrow is facing upward in the insertion direction.</p> <ul style="list-style-type: none"> • Connect the electrode cable for channel I to the appropriate socket (29) on VacoS. • Connect the electrode cable for channel II to the appropriate socket (27) on VacoS. • Insert the red alligator clips to the red connector of the electrode cable. • Insert the black alligator clips to the black connector of the electrode cable. |

7.2 Installation of cables for VacoS

Mounting the Vaco electrode hose.

Connect the black Vaco electrode hose to the appropriate socket (1) on the front of VacoS

at  I - .

Connect the red Vaco electrode hose to the appropriate socket (2) on the front of VacoS

at  I .

Connect the black Vaco electrodes hose to the appropriate socket (3) on the front of VacoS

at  II - .

Connect the red Vaco electrode hose to the appropriate socket (4) on the front of VacoS

at  II + .

Insert the Vaco electrode to the connector on the free end of the Vaco electrode hose.

Assembly Vaco electrodes

Connect each of the 4 Vaco electrode hoses to each Vaco electrode.

Establish network connection

Connect the short mains cable to the sockets provided on VacoS (24) and on Soleo SonoStim / Soleo Galva (14).

Connect the mains cable

Connect the mains cable to the socket provided on VacoS (25) and connect the cable to the mains.

Turn the VacoS on

VacoS is switched on via the toggle switch (23).

Note:

Activation of VacoS take place through this window.



The window is only active when VacoS is on.

8.1 General

Note: *The following descriptions always relate to the treatment with one channel and emanate the factory default settings.*

Note: *Changes to the basic settings can only be made from the start up screen.*

start up Screen Once the device has been switched on and the self-test performed, the start up screen opens.




or



Selecting configuration




Touching the  button opens the "configuration menu".

Configuration menu

In the configuration menu, the factory settings can be changed and adjusted individually. After activating the configuration menu, the "configuration selection" screen becomes active.

8.1 General

Settings General



The settings options are shown below.

By default, the factory settings are pre-programmed as shown in the screen.

Language


Selecting the language



The selection is made directly in the corresponding row.

Program start up settings

Adjustment of the program start up settings.



The selection is made directly in the corresponding row.

Startup screens

Choice of two start up screens.



The selection is made directly in the corresponding row.

8.1 General

Sounds

Switching the signal tones on and off.



The selection is made directly in the corresponding row.

Volume

Adjustment of the volume from 1 to 4.



The adjustment is made using the two arrow keys.

Brightness


Adjustment of screen brightness from 0 to 10.



The adjustment is made using the two arrow keys.

Screensaver

Adjustment of screen saver start time, from 0-20 minutes.



The adjustment is made using the two arrow keys.

Note:

The screensaver function is disabled when therapy is in progress.

Welcome text

Option to enter a personalised welcome text.



Activating the field opens the on-screen keyboard for entering a greeting text.

8.1 General


Set Defaults




Touching button  restores the default factory settings.

Version



Touching the  button opens a window with information regarding the current software version.




Touch the OK button to close the window.

8.2 Electrotherapy

Electrotherapy
settings

Elektrotherapie



Time

Adjustment of a single treatment time of 1 to 60 minutes for all programmes.



The adjustment is made using the two arrow keys.

Polarity reversal cycle Adjustment of periods of time from 10 to 120 seconds.



The adjustment is made using the two arrow keys.

Frequency

Adjustment of the fundamental frequency in stages 2500 Hz, 4000 Hz and 8000 Hz.



The adjustment is made using the two arrow keys.


8.3 Ultrasound Therapy

Note:

Activation of the "SonoSwing" menu and the associated settings are only active in the Soleo SonoStim version.

Ultrasound therapy settings

SonoSwing



Coupling signal


Adjustment (50 to 95%) of the coupling limits value.



The adjustment is made using the two arrow keys.

Unit

Adjustment of unit intensity on the bar graph.



The selection is made directly in the corresponding row.

Maintenance

Wartung

Opens the "Maintenance Configuration Selection" screen.

In the maintenance menu, software updates can be performed.

You will receive the latest information about performing a software update when an update is planned.

To open the maintenance menu, enter the password "armin".

9.1 Electrotherapy

Note:

The following notice applies to all the system therapies.

If the treatment time is extended individually, there may be a change in the method of operation and the patient must be observed during use with increased attention.

Program start up




Touching button  opens the programme screen. Here the desired programme may be selected.

SonoStim programme In SonoStim , 3 different forms of treatment are available:

Selecting electrotherapy

The electrotherapy selection is made directly in the corresponding row.



The following description refers to the electrotherapy programmes.


Selecting the current form group

The selection of the desired current category is carried out directly in the corresponding row (rectangular monophasic currents here).

Soleo SonoStim – Instructions for operation

9

9.1 Electrotherapy




Current form configuration

The current form selection is made directly in the corresponding row. (here NF 00).



Therapy screen


After selection of the current form, the therapy screen opens automatically on channel I.



9.1 Electrotherapy

Therapy startup


By adjusting the intensity from the left intensity regulator, the display changes in the bottom status line from "ready" to "active" and the therapy is being performed. In the bar graph, the current flow is displayed and the therapy time counts down at one-second intervals.




End of therapy

The end of the treatment period is indicated by an acoustic signal at the end of treatment, the clock is set to 00:00, the intensity automatically goes to zero and the bar graph disappears. The display in the lower status bar changes from "active" to "ready".

The treatment time is activated automatically after the end of therapy.




Note:

The vacuum button  is only visible if VacoS is connected to Soleo SonoStim / Soleo Galva as described in section 7.2 and turned on.

9.1 Electrotherapy

Description of the display elements and buttons



(1) Intensity

The bar graph displays the intensity currently selected.

In the CC mode, the display is in mA eff. or peak.

In the CV mode, the display is in peak volts.

The setting is done with the left intensity controller (Mono Stim) or with two controllers (Duostim).

(2) Start / Stop

Pressing the button starts or stops the therapy.

(3) Vacuum

The therapy screen indicates the connection of the vacuum unit.

Activation of the button opens the "vacuum" setting and start-up screen

Note:

If there is no vacuum unit connected, this button is completely hidden.

(4) Channel

After activating the button, several channel mode selection possibilities are available.

In group "mode" you can select between mono/twin and duo.

In the "channel" group, the desired channel for the therapy can be activated.

In the "synchronisation" group, specification of the threshold is determined for the modes twin / duostim. Parallel and separated modes are available alternately.

In duostim mode, the current form can be selected for each channel.

(5) Therapy time

Displays the therapy time for the selected program.

Using the arrow buttons, the preset therapy time can be increased or shortened. After the start of therapy, the current remaining therapy time is displayed.

9.1 Electrotherapy

(6) Status line The status bar shows information regarding the current therapy status. If therapy is not active, the "Ready" display appears; when therapy is in progress the text "Active" appears.

(7) Title bar In the title bar, the name of the current ultrasound therapy program is displayed.

(8) Effect Gives an overview of the medical effects of the current flow form.

Note: If the parameters of the waveform (pulse time or pause) change, the current form's effects will change. Therefore, after a change the window is no longer visible.

(9) Threshold Represents the threshold parameters of the current treatment programme in the activated state. The parameters of increase time, hold time and break time are displayed.
If there is no threshold activated, the threshold window is shown minimised.
Activation of the "threshold" button opens the "threshold parameters" window. Here threshold can be activated or the parameters can be user-defined.
Changed parameters are taken over by activating the "OK" button and the "Cancel" button aborts the process.

Note: If no threshold is provided for a programme, the window is completely hidden.


(10) Parameters Represents the parameters of the current therapy programme.
These are the pulse and pulse pause time (and hence frequency), polarity, CC / CV and galvanic share parameters.
Activation of the "parameters" button opens the "current form parameters" window.
The parameters can be user-defined.
Use the arrow keys to change pulse parameters or the values within the current form up and down.
Change the polarity, change CC/CV and galvanic share by activating the corresponding button.

Note: Modes 1.3 to 1.5 are not useful for all types of current and are therefore not always displayed.

9.2 Ultrasound Therapy


Program start up



Touching button  opens the programme screen. Here the desired programme may be selected.


Selecting ultrasound therapy

The ultrasound therapy selection is made directly in the corresponding row.



Programme selection


The desired ultrasonic therapy program is selected directly in the corresponding row (in this case US 01).



9.2 Ultrasound Therapy

Therapy screen

After selecting the ultrasound therapy program, the therapy screen opens.




Note:


Before starting treatment, check that the data in the parameter window (in this case 5 cm²) corresponds with the connected ultrasound head.


Adjusting the intensity

Adjust the intensity using the left intensity controller.



Therapy startup




Touching  starts the therapy. When therapy starts, the display in the bottom status bar changes from "Ready" to "Active" and the "Start" button changes to "Stop". In the bar graph, the configured dose is displayed and the therapy time counts down at one-second intervals. The coupling display is active.

9.2 Ultrasound Therapy

End of therapy


After expiration of therapy time, an acoustic alert signals the end of therapy and the clock reads 00:00. The intensity automatically returns to zero, the bar graph display disappears and the coupling indicator is inactive. The display in the lower status bar changes from "Active" to "Ready".

The therapy time is automatically re-activated after the end of therapy and the "Stop" button will change to "Start".



9.2 Ultrasound Therapy

Description of the display elements and buttons



(1) Intensity

The bar graph displays the currently set intensity. The intensity is adjusted using the left controller.

(2) Start/Stop

Pressing the button starts or stops the therapy.

(3) Therapy time

Displays the therapy time for the selected program. Using the arrow buttons, the preset therapy time can be increased or shortened. After the start of therapy, the current remaining therapy time is displayed.

(4) Status bar

The status bar shows information regarding the current therapy status. If therapy is not active, the "Ready" display appears. When therapy is in progress the text "Active" appears.

(5) Title bar

In the title bar, the name of the current ultrasound therapy program is displayed.

(6) Coupling

Indicate coupling. The coupling is digitally displayed as a percentage.

(7) Depth effect

The bar graph shows the current frequency ratio, ranging from 0.8 MHz (800 kHz) to 2.4 MHz, as a percentage. The percentage frequency ratio is adjusted using the right controller.


(8) Parameters

Displays the active ultrasound head. Activating the window opens the "Ultrasound parameters" window. Here the ultrasonic mode can be selected and the ultrasound head switched to water bath therapy (see Section 9.3).


9.3 Water Bath Treatment

Note: *If ultrasound therapy is performed in a water bath, the temperature monitoring of the ultrasound head must be modified prior to therapy.*

Performance Activating the parameter button



opens the ultrasound parameters window.



Pressing the water bath button By touching the water bath button



and confirming with "OK", the temperature monitoring of the ultrasound head is modified to suit water bath therapy.

Note: *Upon completion of water bath therapy, it is possible that the temperature of the ultrasound head is too high for therapy outside the water bath. This is displayed in the status bar, with the message "Ultrasound head temperature control". Therapy is not possible during this time. Once the ultrasound head temperature control process is complete, the message disappears and therapy can continue.*

9.4 Simultaneous Therapy


Program start up



Pressing the button **Start** opens the programme screen. Here the desired programme may be selected.


Selecting
simultaneous therapy

The simultaneous therapy selection is made directly in the corresponding row.



Programme selection


The desired simultaneous therapy programme is selected directly in the corresponding row.



9.4 Simultaneous Therapy


Therapy screen

After selecting the simultaneous therapy programme, the therapy screen opens.



Setting the intensity of the ultrasound

Adjust the intensity using the left intensity controller.



Therapy startup

By adjusting the intensity over the right intensity controller, simultaneous therapy is active.

Adjusting electrical stimulation intensity



The display in the lower status bar changes from "ready" to "active". In the left bar graph, the ultrasonic dose appears and the docking indicator is active. The right bar graph shows the actual current flow. The treatment time is reduced by the second.

Note:


In simultaneous therapy it is to be noted that the current-carrying electrode cable (anode) is basically only active for channel I.

9.4 Simultaneous Therapy

End of therapy

After expiration of the treatment time, an acoustic alert signals the end of therapy and the clock reads 00:00. The intensity automatically returns to zero, the bar graph display disappears and the coupling indicator is inactive. The display in the lower status bar changes from "Active" to "Ready".

The treatment time is activated automatically after the end of therapy.




Description of the display elements and buttons

The individual buttons are described in the preceding chapters.

10.1 Electrotherapy

Program start up



Activation of button  opens the electrotherapy programme window.

Selecting the current form group

The selection of the desired current category is carried out directly in the corresponding row (rectangular monophasic currents here).



Note:

The next steps for the implementation of the therapy are described in detail in chapter 9.1.

Soleo SonoStim / Soleo Galva

Instructions for use

11.1 Electrotherapy with VacoS

11


Note:


For the therapy screen, the procedure is described similar to chapter 9.1 electrotherapy.

Note:

If there is no vacuum unit available, the "vacuum" button does not appear. This is also true if the data transfer or the connection of the vacuum unit to Soleo SonoStim/Soleo Galva is disturbed or interrupted.


Opening the vacuum menu


Clicking on the button  opens the vacuum menu.



Activating the basic vacuum

A basic vacuum of 16 kPa is set at the factory.

Pressing the button  activates the basic vacuum for sucking the Vaco electrodes.




With activation of the basic vacuum, the vacuum is active. The bar graph shows the value of the basic vacuum.

Note:

Pulsed mode is set in the factory. Pressing the "mode" button opens the "vacuum" entry field to adjust the pulsed mode.

Selecting the working vacuum



The intensity of the working vacuum is set via the two arrow keys. The bar graph shows the current value of the working vacuum.

Soleo SonoStim / Soleo Galva


Instructions for use

11.1 Electrotherapy with VacoS

11

Starting electrotherapy


By adjusting the intensity from the left intensity regulator, the display changes in the bottom status line from "ready" to "active" and the therapy is being performed. In the bar graph, the current flow is displayed and the therapy time counts down at one-second intervals.



Electrophysiology therapy end

The end of the treatment period is indicated by an acoustic signal at the end of treatment, the clock is set to 00:00, the intensity automatically goes to zero and the bar graph disappears. The display in the lower status bar changes from "active" to "ready".


The treatment time is activated automatically after the end of therapy.



Note:

The vacuum must be manually stopped using "stop" after the time has ended. The "stop" button changes to "start".

Description of the display elements and buttons



Note:

After pressing the "vacuum" button, the parameter window opens.

(1) Start/Stop

Activation or deactivation of this button starts and stops the suction of the electrodes with the basic vacuum.

(2) Mode

Activation of the button opens the "vacuum" entry field for setting the parameters.

(3) Working vacuum

The working vacuum can be adjusted using the arrow keys.

(4) Vacuum display

Indicates the current vacuum in kPa.

(5) Operating mode

Display of the currently selected mode.

(6) Bar graph

Display of the intensity set.

Note:

Pulsed operation is only possible with large vacuum electrodes.

Water separator

VacoS includes an integrated water separator that collects the moisture released from the sponges and thereby protects the vacuum unit from calcification.

If the water separator is full, a message appears on the display "empty water separator".

Note:

If the water separator is not emptied, no therapy can be started.

Empty the water separator as described in chapter 22.

SD card

Custom settings and the list of indications are saved on the SD card.

Note:


If the SD card is not inserted, the message "No SD card found" will appear when you press the "Indications", "Favourites" and "Memory" buttons.

Deactivate the message as described in chapter 22.

The indication menu is used to assist in the choice of therapy.

Selecting a region of the body

A region of the body is selected by clicking the blue circle.



Selecting a symptom

After selecting the desired region of the body (here shoulder), an indications window opens with various disorders in the shoulder area.


The symptom is selected directly in the respective row (in this case shoulder-arm syndrome).



Selecting a differentiated clinical picture


After selecting the symptom, another window opens with differentiated clinical pictures.

The selection of the differentiated disease is carried out directly in the corresponding row (subacute pain here).



Therapy information

After selecting the differentiated clinical picture, another window opens with detailed therapy information and a program recommendation.



Selecting the therapy program



Touching button opens the therapy screen with the corresponding programme.

Soleo SonoStim / Soleo Galva

Instructions for use

11.5 Storing programmes


11


Note:

In order to save programs, the SD card must be properly inserted into the SD card slot. If this is not the case, it is not possible to save programs.

Opening the memory list and naming a program



Touching button  opens the keyboard to enable the entry of the programme name.



There are 2 available options for the programme name:

1. Acceptance of the program name from the input field.
2. Entry of a personal program name. Entry of a personal program name is carried out using the keyboard.

Note:

When entering a personal program name, the program name in the input field must be deleted first.

Note:

The programs can be saved in the favourites list or in the memory list. There are 120 memory slots available.


Soleo SonoStim / Soleo Galva

Instructions for use

11.5 Storing programmes

11

Saving to the
Favourites list /
Memory list



Touching button (1) opens the favourites list and automatically saves the program in the favourites list.

Pressing button (2) opens the memory list and saves the program in the memory list.

Touching the "OK" button closes the "Save" screen and moves the program to the corresponding list.

The program is always stored in the first available slot on the list.

Pressing button (3) cancels the save operation.

Individually saved programs are listed in the favourites list.

Here these can be

1. called up in order to implement therapy or
2. edited (the sequence altered or deleted).

Note:

The steps to retrieve and editing the favourites / memory list are identical; therefore only the retrieve and editing of the favourites list is described.

Selecting the favourites list




Touching button **Favoriten** opens the favourites list.

Calling up a program

The desired program is selected directly in the corresponding row.

Editing the favourites list




By pressing buttons (1) and (2), the individual favourite pages can be viewed. Button (1) scrolls forwards, and button (2) backwards.

Pressing button (3) opens the "Edit Favourites" screen.

Editing favourites

The programme to be processed is selected directly in the corresponding row.



Touching button (1) returns you to the favourites list.

Touching button (2) scrolls the program up.

Touching button (3) scrolls the program down.

Touching button (4) deletes the program.

Note:

Touching button (4) triggers a confirmation prompt:

"Do you really want to delete the program?"


Touching the "Yes" button deletes the program.

Touching the "No" button cancels the delete routine.

Sequence programmes are a combination of up to three electrotherapy programmes that are launched in direct succession one after another automatically. In the programme family "currents for sports rehabilitation", various programmes are predefined at the factory. They usually comprise warming, training or toning and relaxation phases.

Furthermore, there is a freely programmable programme.

Therapy screen
Preset
programme




Note:

Since different pulse shapes and frequencies provide a different current feel, the intensity levels for individual sequences are set in advance of the therapy.

Activating the
sequence window

By activating field this opens the "sequences" window.



Here the intensity thresholds of the individual sequences are defined in advance of sequential therapy.

The sequence selection is made directly in the corresponding row.

Setting intensities

Determining the intensity using the intensity adjuster.

The intensity adjustment is performed according to the sequences 2 and 3.

Note:

*By determining the intensity, the treatment time runs for each sequence.
 By selecting the next sequence the preset intensity of the preceding sequence will be saved automatically.*

With activation of sequence 1, the treatment time becomes active again for the default values and resets the start button.

Therapy startup

Activation of button



starts the therapy.

The 3 sequences take place in succession. The sequence change is indicated by an acoustic signal.

End of therapy

After the treatment period, this is indicated by an acoustic signal at the end of treatment.

12.1 Overview

In Soleo SonoStim / Soleo Galva, there are special programmes for the diagnosis and treatment of paralysis and flaccid paralysis.

With the diagnostic programmes, clues for assessing the severity of flaccid paralysis can be obtained, making it possible in a simple manner to individually set the parameters necessary for the treatment.

Diagnostic programmes

- PA 00 Medium frequency innervation test according to Lange
- PA 01 Faradic test
- PA 02 Chronaxie / accommodation quotient
- PA 15 Gold standard test

For therapy, different programmes are available with the pulse shapes appropriate for paralysis therapy.

The parameters can be individually adjusted according to the disease.


Treatment programmes

- PA 06 Paralysis Treatment with delta currents (alternating)
- PA 07 Paralysis, medium degeneration
- PA 08 Paralysis, low degeneration
- PA 10 Intentional exercises (biphasic)
- PA 11 Muscle rehabilitation force (biphasic)
- PA 12 Muscle rehabilitation endurance (biphasic)

12.2 Medium frequency test by Lange

The medium frequency test according to length gives a rough indication of whether a muscle is partially or completely denervated. The electrodes are placed in the centre of the muscle to be examined.


Display screen



Note:

When using full intensity, the "pulse" button is active.

Performance

Turn on at an appropriate intensity and activating the button  triggers a single pulse.

Subsequently, the test is carried out on healthy muscle. If the two muscles show an approximately equal excitability to the single pulse triggering, the examined muscle is innervated.

If significant differences in the response of muscles to high stimulus intensities are observed, the examined muscle is partially denervated.


12.3 Neofaradic test

The faradic test tells you whether a muscle is innervated normally.

A muscle is normally innervated when it reacts during a stimulation time of 1 minute with an equally long-lasting contraction.

Innervation is to be assumed if no or only short-term contractions are observed.

Display screen




Performance

After selecting the intensity, the stimulus time of 1 minute is started.

12.4 Chronaxy/accommodation quotient

This programme is used to perform a simple and rapid diagnosis. It consists of 3 programme steps in which rheobase, chronaxie and the accommodation thresholds are measured and the accommodation quotient determined.

Display screen



Rheobase

Threshold current to induce a minimal muscle twitch with a rectangular pulse of 1000 ms.

Identifying/assuming rheobase

Step 1: Set to a low intensity
Step 2: Trigger a single pulse

Repeat steps 1 and 2 with increasing intensity until a minimum muscle twitch is visible. With the first minimal muscle twitch, the rheobase value is determined.

Step 3: Accepting the rheobase value

With the acquisition of the detected pulse, the intensity is automatically set to twice the rheobase value. The intensity cannot be changed in this programme step.

Chronaxie

Pulse duration to induce a minimal muscle twitch at a threshold current value, which corresponds to twice the rheobase value.

Identifying/assuming chronaxie

Step 1: Increase (or possibly) decrease the pulse time gradually by activating the arrow keys
Step 2: Trigger a single pulse

Repeat steps 1 and 2 until a minimum muscle twitch is visible. With the first minimal muscle twitch, the chronaxie value is determined.

Step 3: Accepting the chronaxie value

12.4 Chronaxy/accommodation quotient

Accomodation threshold

Threshold current to induce a minimal muscle twitch with a triangular pulse of 1000 ms.

Accommodation quotient

Intensity threshold (mA) of DIC at 1000 ms

Intensity threshold (mA) of RIC at 1000 ms (Rheobase)

**DIC
RIC**

Triangular impulse response
Rectangular impulse response

Identifying/assuming the accomodation threshold


Step 1: Set to a low intensity
Step 2: Trigger a single pulse

Repeat steps 1 and 2 with increasing intensity until a minimum muscle twitch is visible. With the first minimal muscle twitch the accommodation threshold is determined.

Step 3: Accepting the accomodation threshold value

Accommodation quotient

With acceptance of the value of the accommodation threshold, the accommodation quotient is automatically calculated. The measured values are displayed in a window.



12.4 Chronaxy/accommodation quotient

Chronaxie assessment

<0.05 ms hyperexcitability
0.05 to 1 ms normal excitability
1 - 20 ms partial degeneration
> 20 ms full degeneration

Normal chronaxie value (according to Edel)

Anterior muscles (front of the body)
Proximal: 0.08 to 0.16 ms
Distal: 0.16 to 0.32 ms
Posterior muscles (back of the body)
Proximal: 0.16 to 0.32 ms
Distal: 0.44 to 0.72 ms

Assessment of the accommodation quotient according to Edel

6 - 3: normally functioning neuromuscular system
3 - 2: minor/partial degeneration
2 - 1: serious/severe degeneration

12.5 Fisch Gold Test

This programme is used to perform a simple and rapid diagnosis.


The measurement is carried out in 2 phases:

In phase 1, the intensity I_1 of the first muscular contraction is determined at a square pulse of 1 ms pulse width.

In phase 2, the intensity I_2 is determined with the first muscle twitch in a square pulse of 100 ms.

As a result, the quotient $FGQ = I_1 / I_2$ is shown.

Display screen



Performing measurement phase 1

Step 1: Turn on adequate intensity
Step 2: Trigger a single pulse

Repeat steps 1 and 2 with increasing intensity until a minimum muscle twitch is visible.

Step 3: Copy the reading

After the acquisition of the measured value, the parameters are automatically applied to measurement phase 2.

Performing measurement phase 2

Measurement phase 2 is similar to measurement phase 1.

After performing the measurement, the measured values are shown in the top right of the screen.

| | |
|-------------------------|---------------|
| I_1 | 1,5 mA |
| I_2 | 6,2 mA |
| FGQ | 0,25 |

12.5 Fisch Gold Test

Assessment of the quotients

< 1.7: normally functioning neuromuscular system

> 1.7: partial degeneration

Should one of the two intensities not be determined, the degeneration is complete.

13.1 Soleo SonoStim / Soleo Galva

| | |
|-------------------------------|--|
| Power supply | 100 - 240 V / 50 Hz / 60 Hz 220 V / 60 Hz |
| Mains fuse | 2 x T2AL, 250 V, 5x20 mm |
| Electrical energy consumption | max. 60 VA |
| Protection class | I |
| Operating mode | Intermittent operation: 30 minutes on, 10 minutes off |
| Applied Part | Type BF |
| Dimensions | 322mm x 234mm x 130mm |
| Weight | 2.1 kg |
| Operation | +10°C to +35°C, 20 % to 80 % relative humidity without condensation, 700 hPa – 1060 hPa |
| Storage and transport | -10 °C to +50 °C, 10 % to 90 % relative humidity without condensation, at 700 hPa – 1060 hPa |

Note: *Storage and transport in original packaging only.*

Subject to technical modifications!

13.2 Stimulation

Output power

| | | 200 Ω | 500 Ω | 1000 Ω | 2000 Ω |
|--|----|------------------------|------------------------|------------------------|------------------------|
| Galvanisation GA xx | CC | 80 mA | 80 mA | 80 mA | 80 mA |
| | CV | 14 V | 40 V | 60 V | 60 V |
| Diadynamic currents DD xx | CC | 6 mA | 20 mA | 20 mA | 20 mA |
| | CV | 3 V | 8 V | 16 V | 40 V |
| High-voltage currents HV xx | CC | 250 mA _{peak} | 250 mA _{peak} | 250 mA _{peak} | 250 mA _{peak} |
| | CV | 40 V _{peak} | 140 V _{peak} | 200 V _{peak} | 200 V _{peak} |
| Medium frequency currents MF xx | CC | 56 mA | 56 mA | 56 mA | 56 mA |
| | CV | 10 V | 25 V | 60 V | 60 V |
| Interferential currents IF xx | CC | 56 mA | 56 mA | 56 mA | 56 mA |
| | CV | - V | - V | - V | - V |
| Micro stimulation currents MI xx | CC | 1 mA _{peak} | 1 mA _{peak} | 1 mA _{peak} | 1 mA _{peak} |
| | CV | 0.3 V _{peak} | 0.6 V _{peak} | 1 V _{peak} | 2 V _{peak} |
| Low-frequency rectangular currents NF xx | CC | 10 mA _{peak} | 80 mA _{peak} | 80 mA _{peak} | 80 mA _{peak} |
| | CV | 12 V _{peak} | 40 V _{peak} | 80 V _{peak} | 160 V _{peak} |








Polarity of the currents If currents have a polarity: red = positive (+), black = negative (-)

Accuracy $\pm 20\%$

Appliance Stimulation electrodes

13.2 Stimulation

Parameters of current forms

| Current forms | Program mme name | Max. Current (CC) | Frequency | Max. voltage (CV) | Pulse |
|---|------------------------|------------------------|--|-----------------------|--------------|
| Galvanisation GA xx  | GA | 80 mA _{eff} | -- | 60 V _{eff} | -- |
| Diadynamic currents DD xx  | DD | 20 mA _{eff} | DF: 100 Hz, MF: 50 Hz, CP, LP: 50/100 Hz | 60 V _{eff} | 7.5 ms |
| High-voltage currents HV xx  | HV | 250 mA _{peak} | 1-1961 Hz | 128 V _{peak} | 10 us-0.1 ms |
| Medium frequency currents MF xx  | MF | 80 mA _{peak} | 8000/3-200 Hz | 80 V _{peak} | 1 ms-50 ms |
| Interferential currents IF xx  | IF | 100 mA _{peak} | 8000/10-120 Hz | -- | 125 us-0.4 s |
| Micro stimulation currents MI xx  | MI | 2 mA _{peak} | 2Hz, 100 Hz | 20 V _{peak} | 20 us-1s |
| Low-frequency rectangular currents NF xx  | NF | 80 mA _{peak} | 1-1923 Hz | 180 V _{peak} | 20 us-1s |

All listed parameters are valid for a load resistance of 1 KOhm.

(Pulse duration: 20 us - 1 s, pulse repetition frequency: 1-1923 Hz, pulse repetition amplitude: 80 mA_{peak}, DC voltage: 50 V. These values are for a load resistance of 500 ohm)

13.3 Ultrasound

Ultrasound heads

| | |
|-------------------------------|---|
| Frequency | 0.8 MHz and 2.4 MHz |
| Ultrasound head, small | 1 cm ² , ERA = 0,67cm ² at 0.8 MHz, 0.65 cm ² at 2.4 MHz |
| Maximum power | 1.0 W at 0.8 MHz, 0.6 W at 2.4 MHz |
| Intensity levels | 0.1 to 1 W/cm ² eff. at intervals of 0.1 W / cm ² |
| Ultrasound head, large | 5 cm ² , ERA = 2.30 cm ² at 0.8 MHz, 2.38 cm ² at 2.4 MHz |
| Maximum power | 6.9 W at 0.8 MHz, 7.1 W at 2.4 MHz |
| Intensity levels | 0.1 to 3 W/cm ² eff. at intervals of 0.1 W/cm ² |
| Accuracy | ± 20 % |
| Sound forms | 1. Continuous ultrasound 2. Impulse ultrasound, adjustable pulse frequencies: 20 Hz, 50 Hz, 100 Hz Duty ratio: 1: 1, 1: 2, 1: 3, 1: 5, 1: 10 |
| Interchangeability | Ultrasound heads are calibrated ex works and can be easily replaced. |
| Applied Dart | Transducer ultrasound head |

13.4 VacoS

| | |
|-------------------------------|---|
| Mains voltage | 100 – 120 V / 50 Hz / 60 Hz (1) 230 – 240 V / 50 Hz / 60 Hz 220 V / 60 Hz (2) |
| Mains fuse | 2 x T 0.4 AL, 100 – 120 V, 5 x 20 mm 2 x T 0.2 AL, 230 – 240 V, 5 x 20 mm |
| Electrical energy consumption | max. 50 VA |
| Protection class | I |
| Applied Part | Type BF |
| Operating mode | Intermittent operation: 30 minutes on, 10 minutes off |
| Dimensions | 322mm x 234mm x 130mm |
| Weight | 3.0 kg |
| Operation | +10°C to +35°C, 20% to 80% relative humidity without condensation, 700 hPa – 1060 hPa |
| Storage and transport | -10 °C to +50 °C, 10% to 90% relative humidity without condensation, at 700 hPa - 1060 hPa |
| Negative pressure | 12-60 kPa |
| Pulse operation | Period adjustable to 1 to 8 seconds in steps of 0.5 seconds, ratio 1: 1 |
| Accuracy | ± 20% |
| Applied Part | Vacuum suction cups and sponges |

Note: *Storage and transport in original packaging only.*

Subject to technical modifications!



- Before starting any cleaning and maintenance measures, the device must always be switched off at the main switch and unplugged.
- Make sure that during cleaning and disinfection no liquids penetrate the device. Do not use sprays.
- If liquid penetrates the device during cleaning or disinfecting, please put the unit out of service, protect it from getting used again and contact your service representative.
- Make sure that when cleaning and disinfecting the labels of the device (such as warning, labels of control devices, indication plate) are not damaged.
- The device and its applied part are not considered critical in relation to hygiene when used on non-injured and healthy skin.

Housing

Cleaning: The housing and all cables can be cleaned using commercially available, alcohol-free cleaning agents designed for plastics. Wipe the surface until the dirt is removed using a soft cloth soaked according to the specifications of the manufacturer of the cleaning agent so it is not dripping.

Disinfection: We recommend disinfecting the device at least once a week or upon signs of contamination. Consult with your health professional when doing so. Always perform cleaning prior to disinfection.

Housing can be disinfected using disinfectant wipes. To do so, use a commercially available alcohol-free disinfectant for metal and plastic, with bactericidal, virucidal and fungicidal properties. Observe the application instructions of the manufacturer. Wipe all surfaces using a soft cloth soaked according to the specifications of the manufacturer of the disinfectant, but not dripping, or use pre-impregnated disinfectant wipes (so called wipes). If applicable, also observe requirements for drying or post-cleaning.

Ultrasound heads

Cleaning: Proceed as described under "housing".

Disinfection: Proceed as described under "housing".

Vaco electrodes

Cleaning: The Vaco electrodes can be cleaned of any visible dirt with commercially available soft plastic cleaners. Unplug the Vaco electrode from the electrode hose and invert the suction cup completely. Dip the Vaco electrode into a cleaning solution prepared according to the manufacturer of the detergent and clean the electrode thoroughly using a stiff brush. Make sure that the gap underneath the metal plate is cleaned completely safely. Rinse with clear water.

Disinfection: We recommend disinfecting the device at least once a week or upon signs of contamination. Consult with your health professional when doing so. Always perform cleaning prior to disinfection.

The Vaco-electrodes can be disinfected with a commercial alcoholic disinfectant for metal and plastic, which is bactericidal, virucidal and fungicidal. Observe the application instructions of the manufacturer. Dip the Vaco electrodes fully in the disinfectant, and move around a little. Make sure that the inner and outer surface is completely covered with disinfectant.

If applicable, also observe requirements for drying or post-cleaning.

Sponges/sponge bags

Cleaning: Rinse the sponge/sponge bags thoroughly with water. The use of cleaning agents is not recommended as components remain in the sponges and can cause skin irritation during application. Furthermore, it may lead to damage of the equipment.

Disinfection: For thermal disinfection, heat the sponges/sponge bags for 10 minutes or wash (water temperature > 95°C) in pure water without additives. The use of disinfectants is not recommended because the elements that remain in the sponges may result in skin irritation during application. In addition, the equipment may be damaged.



If, however, flammable solutions are used for cleaning and disinfection, enough time must be allowed for them to evaporate before using the device. This may otherwise result in ignition.

Note:

Only use the device in a perfectly hygienic environment.

General note:

Selection and application of the electrodes must be done with care. In constant-current operation, ensure good and uniform contact with the skin. A reduction in the contact area can lead to discomfort for the patient.

Before treatment, the skin should be inspected and cleaned if necessary, for example if the patient is sweaty or ointments have been applied. Inflamed skin, small wounds or cracks are to be covered with Vaseline or zinc paste. Exercise extreme caution with recent scars.

Single-use electrodes

The self-adhesive disposable electrode allows for convenient and rapid application. Through one-time use, it is hygienic to use. Three therapeutically useful sizes allow individual therapy tailored to the clinical picture.

Single-use electrodes are particularly suitable for treatment with bipolar flows. For therapy with monopolar pulse currents or currents with a galvanic component, the disposable electrodes should be additionally padded with a damp sponge.

For the Soleo SonoStim/Soleo Galva, all electrodes standard in electrotherapy must be connected by means of an insulated cable clamp.

Note:

The disposable electrode is reserved for one-time use and can easily be disposed of as household waste.

Repeated use of the disposable electrode can be dangerous for the patient.

For particularly large-scale applications in addition to the disposable electrodes, plate and rubber electrodes are also available.

The electrical connection is made as in the disposable electrode to the terminals of the electrode cables. Attach clips to the plate electrode fully inserted into the sponge bag.

Rubber electrodes

Rubber electrodes are suitable for therapy with bipolar currents. In the application of pure galvanic current flows with electrical component or long pulse duration, it is important to note that through normal use conditional removal of carbon results in a reduction of the conductivity.

Tin plate electrodes

For the treatment with pure galvanic current (galvanisation, iontophoresis) large tin plate electrodes are suitable.

Sponge bags and sponges

Both rubber and tin plate electrodes must always be backed by a wet liner. For rubber electrodes, the use of sponge bags is recommended for tin plate electrodes, which should be at least 1 to 2 cm thick. To moisten the sponges, tap water is recommended. Distilled water is not suitable, because of its poor conductivity.

In contrast to the convenient self-adhesive properties of single-use electrodes, rubber and tin plate electrodes must be fixed. For this purpose, Velcro or rubber bands are suitable.

Fully insert the electrodes into the sponge bags and fit using slight pressure until they rest properly on the body. The straps should not leave any constrictions.

The products bear the CE marking



in accordance with EC Directive 93/42/EEC concerning medical devices.

Manufacturer

Zimmer MedizinSystems GmbH
Junkersstraße 9
89231 Neu-Ulm, Germany
Tel. +49 731. 9761-291
Fax: +49 731. 9761-299
www.zimmer.de

Scope of delivery
item no.

5302
(see below)

Soleo SonoStim

1 basic device Soleo SonoStim
2 pairs of electrode cables with twist protection and magnetic clip,
incl. 2 red and 2 black alligator clips, 2.90 m long
1 ultrasound head 0.8 and 2.4 MHz, ø 28 mm
1 tray mount, right
1 tray mount, left
1 mains cable*
2 touch pens
2 test resistors
1 instructions for use

Accessories
item no.

65800410 Touch pen
118 Mains cable*
65910321 Tray mount, right
65910311 Tray mount, left
87054001 Pocket for Soleoline
87053000 Case for Soleoline
10101393 Instructions for use

Electrotherapy

154 1 pair of electrode cable with twist protection and magnetic clip,
incl. 1 red and 1 black alligator clip, 2.90 m long
68910911 1 pair of electrode cables with twist protection and magnetic clip, 2.90 m long
31100146 Alligator clip, red
31100147 Alligator clip, black
87200120 Box single-use electrodes, small, 32 x 40 mm, 200 pairs
87200140 Box single-use electrodes, medium, 56 x 56 mm, 100 pairs
87200130 Box single-use electrodes, large, 56 x 128 mm, 50 pairs
87200110 Box single-use electrodes, round, ø 100 pairs
43 Electrode plate of pure tin, 90 x 120 mm
96 Sponge for electrode plate 120 x 160 mm
212 1000 films for iontophoresis, 130 x 180 mm
44 Rubber electrode 50 x 50 mm, 1 pair
97 Sponge pocket 70 x 75 mm for # 44, 1 pair
46 Rubber electrode 50 x 100 mm, 1 pair
98 Sponge pocket 115 x 80 mm for # 46, 1 pair
232 Rubber strap 60 cm, with button
233 Rubber strap 120 cm, with button
230 Velcro tape, 60 cm long, 10 cm wide
231 Velcro tape, 120 cm long, 10 cm wide
65920310 Test resistor

Ultrasound therapy

Art. No.
4200 Ultrasound head 0.8 and 2.4 MHz, ø 28 mm
4220 Ultrasound head 0.8 and 2.4 MHz, ø 13 mm
6 Sono plus, 1 bottle

Scope of delivery

Art. No.

5330

(see below)

Soleo Galva

- 1 basic device Soleo Galva
- 2 pairs of electrode cables with twist protection and magnetic clip,
incl. 2 red and 2 black alligator clips, 2.90 m long
- 1 tray mount, right
- 1 tray mount, left
- 1 mains cable*
- 2 touch pens
- 2 test resistors
- 1 instructions for use

Accessories

Art. No.

65800410

118

65910321

65910311

10101393

- Touch pen
- Mains cable*
- Tray mount, right
- Tray mount, left
- Instructions for use

154

- 1 pair of electrode cable with twist protection and magnetic clip,
incl. 1 red and 1 black alligator clip, 2.90 m long

68910911

- 1 pair of electrode cables with twist protection and magnetic clip, 2.90 m long

31100146

- Alligator clip, red

31100147

- Alligator clip, black

87200120

- Box single-use electrodes, small, 32 x 40 mm, 200 pairs

87200140

- Box single-use electrodes, medium, 56 x 56 mm, 100 pairs

87200130

- Box single-use electrodes, large, 56 x 128 mm, 50 pairs

87200110

- Box single-use electrodes, round, ø 100 pairs

43

- Electrode plate of pure tin, 90 x 120 mm

96

- Sponge for electrode plate 120 x 160 mm

212

- 1000 films for iontophoresis, 130 x 180 mm

44

- Rubber electrode 50 x 50 mm, 1 pair

97

- Sponge pocket 70 x 75 mm for # 44, 1 pair

46

- Rubber electrode 50 x 100 mm, 1 pair

98

- Sponge pocket 115 x 80 mm for # 46, 1 pair

232

- Rubber strap 60 cm, with button

233

- Rubber strap 120 cm, with button

230

- Velcro tape, 60 cm long, 10 cm wide

231

- Velcro tape, 120 cm long, 10 cm wide

65920310

- Test resistor

* Standard cable. Other country-specific plug types are available. If necessary, contact your dealer.

Scope of delivery

Accessories

17.3 VacoS

17

Scope of delivery
Art. No.
5320

VacoS

1 basic device VacoS
1 pair Vaco electrode tubes red, 2.10 m
1 pair Vaco electrode tubes black, 2.10 m
2 pairs Vaco electrodes ø 90 mm, self-closing with sponges
2 pairs Vaco electrodes ø 60 mm, self-closing with sponges
2 connection cables for Soleo *SonoStim* / Soleo *Galva*
1 mains connection cable, 50 cm
4 pieces sorting combs for Vaco lines incl. magnetic clip
1 water bottle

Accessories
Art. No.

| | |
|----------|---|
| 95 | Sponge for Vaco electrode, large, for # 72 |
| 94 | Sponge for Vaco electrode, small, for # 71 |
| 72 | Vaco electrode, ø 90 mm, self-closing |
| 71 | Vaco electrode, ø 60 mm, self-closing |
| 164 | Vaco electrode tube, red, 2.10 m |
| 165 | Vaco electrode tube, black, 2.10 m |
| 68910210 | Connection cable for Soleo <i>SonoStim</i> / Soleo <i>Galva</i> |
| 119 | Mains connection cable, 50 cm |
| 65350110 | Sorting comb for Vaco lines incl. magnetic clip |
| 95920000 | Water bottle |

Soleo *SonoStim* and Soleo *Galva* can be used with *VacoS*.

If you combine the devices and therefore operate a medical system, you are solely responsible for the accuracy of the connection.

Soleo *SonoStim* / Soleo *Galva* are manufactured according to EN 60601-1 safety regulations.

As the manufacturer, Zimmer MedizinSysteme can only be considered responsible for safety and reliability if

- the equipment is operated using a proper power outlet with an earthing contact and the electrical installation complies with DIN VDE 0100 Part 710,
- the device is operated in accordance with the user manual,
- extensions, readjustments or modifications are only carried out by persons authorised by Zimmer MedizinSysteme,
- the user is satisfied regarding the functional safety and the proper condition prior to using the device,
- the ultrasound head, cables and connectors are inspected for any damage that may affect the safety of the device (e.g. cracks) prior to use,
- the equipment is only operated by properly trained personnel,
- the device is immediately disconnected from the mains in the event of the penetration of liquids,
- the device is not operated in hazardous areas and / or a combustible atmosphere.

The device does not contain any parts that can be serviced and repaired by the operator.



Fuses and other spare parts may be only be exchanged by trained service personnel. Replacement of the lithium battery by untrained persons may cause risks.



Servicing of the device must be performed by trained personnel. All information necessary for service is provided in the Soleoline service manual or can be obtained from the manufacturer. On request, Zimmer MedizinSysteme provides circuit diagrams, component lists, descriptions, calibration instructions or other documents.

Soleo *SonoStim* / Soleo *Galva* performs a self-test after turning on, which checks all internal components.

In the event of a fault, an error message appears.

In addition, an extended function test can be performed as described below for all 3 modes.

These tests should be performed monthly or in case of any doubt regarding the correct function of the device.

Electrical muscle stimulation

Selecting programme GA 00. Connecting the test adapter.
At maximum intensity, the display in the bar graph must be full. Run the test with both channels, one after the other.

Ultrasound

Select the ultrasound head and cover the ultrasound head surface with coupling gel. In the case of starting therapy with low power, the coupling display must show more than 90%.
Run the test with both ultrasound heads, one after the other.
Then clean the ultrasound heads.

Vacuum

Select vacuum and create a basic vacuum. Short circuit the Vaco electrode and set the maximum vacuum. In non-pulsed operation, the bar graph must be full.

Vaco electrodes

Vaco electrodes evert outwardly so that the electrode area is no longer enclosed by the rubber-lip and a vacuum is created.
VacoS must suction permanently.

The Soleo *Galva*, Sono*Stim* and optionally available Vaco*S* devices are listed in Annex 1 of MPBetreibV (Medical Device Operating Directive). Please observe the measures stipulated therein.

The devices are not listed in Annex 2 of MPBetreibV (Medical Device Operating Directive).

In Germany, the DGUV regulation 3 (Electrical Systems and Equipment) in its current version must also be observed.

Note:

These notes apply to the operation of the device in Germany. If applicable, observe any divergent national regulations that apply in your country.

Cable inspection



Means with the current stimulation therapy in constant current mode:
Interruption of the patient circuit.

As a rule, this message concerns disconnected electrodes, dirty electrode terminals, a defective patient cable or defective patient safety. Clear message by confirming "OK".

Overcurrent



This refers to an increase in the maximum allowable current. An indicated increase in current in constant operation usually indicates a fault in the device, while increasing the current in constant voltage mode by changing the patient's resistance may occur (e.g. wet skin).

Clear message by confirming "OK".

If the error message occurs again, please contact customer service.

No SD card found



If the SD card is not inserted, the message "No SD card found" will appear when you press the "Indications", "Favourites" and "Memory" buttons.

Insert card and confirm with "OK".

VacoS

No vacuum

If the *VacoS* *should* not reach the set vacuum, there may be a leak in the system.

Check the following in this case,

- whether all four Vaco connected electrode hoses are closed and connected to a Vaco electrode
- whether the Vaco electrodes are applied to the patient properly.

Vacuum button disabled

Should the vacuum button be disabled, a transmission error may be present.

Check the following in this case,

- whether Soleo *SonoStim* / Soleo *Galva* is connected correctly to *VacoS*
- whether Soleo *SonoStim* / Soleo *Galva* is connected to *VacoS* via the connecting cable
- whether *VacoS* is turned on.

Note:

VacoS and Duostim mode isolated

In "isolated Duostim" mode, a double-circuit vacuum application is not possible. In "Duostim mode", only channel I is active for the vacuum application. Channel II cannot be operated with vacuum application.

VacoS

Water separator full

If the water separator is full, a message appears on the display "empty water separator".
The water is emptied by means of a water bottle that comes with VacoS. Close the message "empty water" by pressing "OK".

Note:

*If the water separator is not emptied, no therapy can be started.
We recommend emptying the water separator daily.*

Empty the water separator

Mounting the bottle

The connection valve (33) for draining the water separator is on the right side of VacoS. On the water bottle hose, there is a quick-release lock which is connected to the connection valve on the device.

Empty the water

Squeezing the bottle creates a vacuum, which empties the water separator. Under certain circumstances, a repeated structure of the negative pressure is necessary to empty the water separator.


Dismantling/emptying the bottle

By squeezing the sides of the quick connector on the tube of the water bottle and simultaneously removing the terminal, the bottle is removed from the valve connection.

The lid of the bottle is removed by turning so that the bottle can be emptied.

Error messages

A message regarding the error that has occurred appears on the display, e.g.:



Some errors are also displayed with a corresponding error number (e.g. "Self test error X/XX/X").

In some cases, the error can be corrected by switching off, waiting for five seconds and then switching the device on again.

Otherwise, please contact customer service.

You can reach the service via your assigned sales representative or our headquarters in Neu-Ulm.

The device must be returned to the factory in its original packaging.

For questions or equipment malfunctions, please contact the address listed:

Headquarters

Zimmer MedizinSysteme GmbH
Junkersstraße 9
89231 Neu-Ulm, Germany
Tel. +49 731. 9761-291
Fax: +49 731. 9761-299
www.zimmer.de

Disposal

The device may only be returned to the factory in its original packaging.
It must only be disposed of via the factory in Neu-Ulm.

In foreign (European) countries please refer to national regulations for disposal.
Contact your distributor if necessary.

Medical electrical equipment such as the Soleo SonoStim/Soleo Galva are subject to special precautions with regard to EMC (electromagnetic compatibility) and must be installed and commissioned in accordance with the EMC information contained in the user manual or the accompanying documents.

Portable and mobile RF communications equipment (e.g. mobile phones) can affect electrical medical equipment.

Soleo SonoStim / Soleo Galva may only be operated with the original parts specified in the list of the delivery and accessories. The operation of the device with other parts may result in increased emissions or in the device's decreased interference immunity.

| Guidelines and Manufacturer's Declaration - Electromagnetic Emissions | | |
|--|------------|---|
| The Soleo SonoStim / Soleo Galva device is designed for operation in the electromagnetic environment specified below. The customer or the user of the Soleo SonoStim/Soleo Galva must ensure that it is operated in such an environment. | | |
| Emissions measurements | Compliance | Electromagnetic environment - guidelines |
| RF emissions in accordance with CISPR 11 | Group 2 | The Soleo SonoStim/Soleo Galva device must emit electromagnetic energy in order to ensure its intended function. Nearby electronic equipment may be affected. |
| RF emissions in accordance with CISPR 11 | Class B | |
| Emissions of harmonics in accordance with IEC 61000-3-2 | Category A | |
| Emissions of voltage fluctuations / flickers in accordance with IEC 61000-3-3 | Compliant | The Soleo SonoStim/Soneo Galva device is suitable for use in all establishments, including domestic establishments and those directly connected to the public supply network that also supplies buildings used for domestic purposes. |


Table 201 according to EN 60601-1-2: 2006-10

The device must not be used in direct proximity to or stacked directly on top of another similar device. If operation near to or stacked on top of another device is unavoidable, the device should be monitored to verify its proper operation within this setup.

| Guidelines and manufacturer's declaration - electromagnetic immunity | | | |
|---|--|--|--|
| The Soleo SonoStim/Soleo Galva device is designed for operation in the electromagnetic environment specified below. The customer or the user of the Soleo SonoStim/Soneo Galva device should ensure that it is used in such an environment. | | | |
| Immunity tests | IEC 60601 - Test Level | Compliance level | Electromagnetic environment - Guidelines |
| Electrostatic discharge (ESD) in accordance with IEC 61000-4-2 | ± 6 kV contact discharge ± 8 kV air discharge | ± 6 kV contact discharge ± 8 kV air discharge | Floors should be made from wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity must be at least 30%. |
| Electrical fast transient / burst in accordance with IEC 61000-4-4 | ± 2 kV for mains cables ± 1 kV for input and output cables | ± 2 kV for mains cables Not applicable | The supply voltage quality must correspond to that of a typical commercial or hospital environment. |
| Surges in accordance with IEC 6100-4-5 | ± 1 kV differential mode ± 2 kV common mode | ± 1 kV differential mode ± 2 kV common mode | The supply voltage quality must correspond to that of a typical commercial or hospital environment. |
| Voltage dips, short-term interruptions and voltage variations in accordance with | < 5% U_T (> 95% dip in U_T for ½ period) 40% U_T (60% dip in U_T for 5 periods) 70% U_T (30% dip in U_T for 25 periods) < 5% U_T (> 95% dip in U_T for 5 seconds) | < 5% U_T (> 95% dip in U_T for ½ period) 40% U_T (60% dip in U_T for 5 periods) 70% U_T (30% dip in U_T for 25 periods) < 5% U_T (> 95% dip in U_T for 5 seconds) | The supply voltage quality must correspond to that of a typical commercial or hospital environment. If the user of the Soleo SonoStim/Soneo Galva requires continued operation, even in the case of interruptions in the power supply, it is recommended that the Soleo SonoStim/Soneo Galva be powered from an uninterruptible power supply or a battery. |
| Magnetic field of the supply frequency (50/60 Hz) in accordance with IEC 61000-4-8 | 3 A/m | 3 A/m | Magnetic fields at mains frequency should have the typical values found in a business or hospital environment. |
| Note: U_T is the mains AC voltage before application of the test level. | | | |

Table 202 according to EN 60601-1-2: 2006-10

Main features of the Soleo SonoStim / Soleo Galva (only electrical stimulation) are: trouble-free delivery of electrical stimulation and ultrasound, in conjunction with the vacuum unit and vacuum with the set parameters and trouble-free operation of all functions.

| Guidelines and manufacturer's declaration - electromagnetic immunity | | | |
|---|---|--|--|
| The Soleo SonoStim/Soleo Galva device is designed for operation in the electromagnetic environment specified below. The customer or the user of the Soleo SonoStim/Soneo Galva device should ensure that it is used in such an environment. | | | |
| Immunity tests | IEC 60601 – Test Level | Compliance level | Electromagnetic environment - Guidelines |
| <p>Conducted RF disturbance variables in accordance with IEC 61000-4-6</p> <p>Radiated RF disturbance variables in accordance with IEC 61000-4-3</p> | <p>3 V_{RMS} 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p> | <p>3 V_{RMS} 150 kHz to 80 MHz</p> <p>10 V/m 80 MHz to 2.5 GHz</p> | <p>Portable and mobile RF communications equipment should not be used any closer to the Soleo SonoStim/Soneo Galva, including cables, than the recommended separation distance calculated using the equation relevant for the transmission frequency.</p> <p>Recommended separation distance:</p> <p>$d = 1.17 \sqrt{P}$</p> <p>$d = 0.35 \sqrt{P}$ for 80MHz to 800 MHz</p> <p>$d = 0.7 \sqrt{P}$ for 800 MHz to 2.5 GHz</p> <p>where P is the power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>In accordance with an on-site investigation^a, for all frequencies, the field strength of a stationary radio transmitter should be less than the coincidence level^b.</p> <p>Errors may occur in the vicinity of devices marked with the following symbol:</p>  |
| <p>NOTE 1: At 80 Hz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | |

^aThe field strengths of fixed transmitters e.g. base stations of mobile telephones and mobile radio devices, amateur radio, AM and FM radio broadcasting, and TV broadcasting cannot be accurately predicted. In order to assess the electromagnetic environment with regard to fixed RF transmitters, a study of the electromagnetic phenomena at the location should be considered. If the measured field strength in the location in which the Soleo SonoStim/Soleo Galva device is used exceeds the above compliance level, the Soleo SonoStim/Soleo Galva device must be monitored to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as the realignment or relocation of the Soleo SonoStim/Soleo Galva device.

^bIn the frequency range of 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Soleo SonoStim/Soleo Galva device

The Soleo SonoStim/Soleo Galva device is intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or the user of the Soleo SonoStim/Soleo Galva device can thus help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Soleo SonoStim/Soleo Galva device – depending on the output power of the communications equipment as indicated below.

| Rated power of the transmitter W | Separation distance depending on the frequency of transmitter m | | |
|-------------------------------------|--|-----------------------------------|--|
| | 150 KHz to 80 MHz $d = 1.17 \sqrt{P}$ | 80–800 MHz $d = 0.35 \sqrt{P}$ | 800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$ |
| 0.01 | 0.12 | 0.04 | 0.07 |
| 0.1 | 0.37 | 0.11 | 0.22 |
| 1 | 1.17 | 0.35 | 0.7 |
| 10 | 3.70 | 1.11 | 2.21 |
| 100 | 11.67 | 3.5 | 7.0 |

For transmitters for which the maximum nominal output is not specified in the above table, the recommended separation distance d in metres (m) can be determined using the equation given for the corresponding column, where P is the maximum rated power of the transmitter in watts (W) in accordance with the manufacturer of the transmitter.

NOTE 1: At 80 and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Soleo*line*

Instructions for Use

Zimmer MedizinSysteme GmbH
Junkersstraße 9
D-89231 Neu-Ulm
Tel. +49 7 31. 97 61-291
Fax +49 7 31. 97 61-299
export@zimmer.de
www.zimmer.de

Zimmer
MedizinSysteme