




PHYSIOTHERAPY & MASSAGE RILA

Ref. EL032.1.A66.SH.WF

5-section physiotherapy and osteopathy couch for a wide range of treatments. Equipped with armrests on the head section, adjustable in height to ensure optimal positioning when the patient is lying in the prone position. Headrest adjustable in positive and negative angles for cervical and cranial manipulation. Fixed middle section. Roll holder included as standard.



Watch video



Specifications:

- Electric height adjustment via perimeter control.
- Elegant design with rounded feet for enhanced safety.
- Headrest with adjustable armrests.
- Headrest with face hole and plug.
- Levelling feet for uneven surfaces.
- Retractable castors ø 75mm.
- High-reliability TiMOTION motor.
- Lifting time: 30 seconds.
- Padding thickness: 5 cm.
- Foam density: 36 kg/m³
- Fire-retardant, highly wear-resistant vinyl upholstery.
- Certifications: CE, FDA, ISO 13485:2016
- Load capacity: 250 kg
- Upholstery colours (options: Blue or Grey)
- Warranty: 1 year for frame, electrical components and upholstery (except misuse).



65.0



0.72



194.0x70.0x53.0 cm