

Fiji - Pacific examination and treatment tables



Product short description:

Fiji - examination and treatment table is a modern 5 and 7-section table that features **electric height adjustment, a heated table surface and LED lighting - (optional)**. It is reliable, aesthetically pleasing, and comfortable to use. The table's modern design will enhance the prestige of your practice and make your clients feel valued.

Product gallery:





Product description:

Fiji - Pacific examination and treatment table

This modern five and seven section examination and treatment table features an ergonomic design that ensures therapist comfort. The wide lying surface ensures patient safety and comfort.

Features:

Five or seven section examination and treatment table with electric height adjustment by the means of a remote (controller).

Head section adjustment by the means of gas spring.

Head rest with shaped face hole with crescent plug.

Foldable armrests lowered till 19 cm (only in 7 section table).

Adjustment of: back, thighs, shanks sections by means of electric actuators.

Multi-layered and very soft padding of 15 cm thickness with upholstery in many colours.

Trendelenburg/reverse Trendelenburg position (15°).

Two castors

Additional options:

An easy-to-use heating panel with two temperature settings for the padding of the **Fiji examination and treatment table**.

Trendelenburg/reverse Trendelenburg position (15°).

Two castors

The remote control is designed to be simple and easy to use:

Fiji - examination and treatment table is equipped with a 10-button remote control that allows quick repositioning of the chair:

- Backrest adjustment
- Thigh section adjustment
- Shank section adjustment
- Height adjustment of the table
- Tilt table (chair function)

Fiji is a luxury examination and treatment table covered with unique upholstery that is:

- Stretch-resistant
- Extremely soft and velvety
- Easy to clean
- Covered with **Silverguard®** for fungal and bacterial resistance.
- Resists many liquids, including **alcohol, oil, blood, and urine**.
- Low temperature resistant
- UV resistant
- Abrasion resistant

The remote control is designed to be simple and easy to use:

Numebr of sections	Height adjustment [cm]:	5 or 7 62-
Height adjustment of armrests [cm]:	Headrest	10259-99 0-
section angle adjustment [°]:	Chest section	19
angle adjustment [°]:	Thigh section angle	0/-500/+38
adjustment [°]:	Shank section angle adjustment	0/+42 0/+25
[°]:	Trendelenburg / anti-Trendelenburg position	0/+18
[°]:	Upholstery thickness [cm]:	0/+150/-15 9
		or 15

Dimensions (L x W) [cm]:	212 x 80 150
Weight [kg]:	150
Max. lifting capacity [kg]:	220-240 / 50-60
Power supply [V/Hz]:	

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Type of table: Examination table
 Type of table: Treatment table
 Number of sections: 5 Number of sections: 7
 Wheeling system: F0 Length [cm]: 212 Width [cm]: 80
 Weight [kg]: 150 Max. lifting capacity [kg]: 150
 Power supply [V / Hz]: 230 / 60 Warranty: 24 months
 Medical device: Yes

Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange
 Tables - Fiji / Hawaii: 5 - section treatment table, 7 - section treatment table
 Tables - Additional options: -, Heated table surface + LED Lightning