

Theraspace Z3 - Massage table

Product short description:



Theraspace Z3 is a professional massage table, designed for demanding physiotherapists. Its modern design and functionality make it an ideal tool for everyday work.

The upholstery on the table is covered by a 36-month warranty.

See more details on product page.

Product gallery:





Product description:

Theraspace Z3 - Massage and rehabilitation table

The Theraspace massage table is designed for performing manual massage treatments to treat and alleviate traumas or disabilities by restoring or improving the function of massaged tissues, organs and systems. The clinical benefits of massage tables are based on the fact that it is possible for the patient to take the right position for the massage treatment, the possibility to adjust the height and/or individual sections of the table, which facilitates the work of the rehabilitator/physiotherapist and increases the comfort of the patient during the massage treatment.

Theraspace Z3 - massage and rehabilitation table - appreciated by therapists all over the world.

Features:

3-section massage table.

High-quality, soft and pleasant to the touch polyurethane upholstery.

Progressive foam.

Electric adjustment of: head + backrest, seat + leg sections via control switches.

Two-layer padding (3 cm foam 90 kg/cm³ + 1 cm foam 43 kg/cm³).

Round breathing hole in head section with 2 types of plugs (crescent, round).

High/low height adjustment foot bar (49-95 cm).

Total length of the bed 205 cm.

2 traction support sleeves

Protection against accidental start-up.

Linak actuator.

LINAK® Technology

Theraspace tables are based on LINAK® components using the latest actuators with IRS technology. This results in increased safety and minimised risk of damage to the equipment.

SmartS4Arms

SmartS4Arms is a system for smooth, screwless adjustment of the armrests. The height is adjusted by palm pressure, and an automatic position lock prevents unwanted movement of the armrests during treatment.

Theraspace Series

A unique experience for your patients

For the patient, comfort during the treatment is as important as the results. With the Theraspace Series of massage tables, you can guarantee your patients comfort and a sense of security during the therapy!

The Theraspace Series massage table features:

Two-layer padding (3 cm foam 90 kg/cm³ + 1 cm foam 43 kg/cm³).

Pivot option.

Stretched upholstery.

A frame combining wood and steel.

F4 wheel system consisting of 4 liftable directional wheels with a central locking mechanism and 4 rubber feet.

High-quality, soft upholstery in a wide variety of colors.

Theraspace - Massage tables are widely used in:

Sanatorium

Massage tables from the Theraspace Series are a top-quality choice for sanatoriums, offering a comprehensive range of solutions designed to meet the needs of even the most demanding customers. Built with advanced technology and precision craftsmanship, these tables fully **support the recovery and relaxation process** of the client.

Rehabilitation clinics

Theraspace massage tables are unparalleled in their class, providing comprehensive solutions for rehabilitation centres that are critical to the recovery process and **improving patients' quality of life**. These tables offer not only comfort but also **effective therapy**, thanks to their advanced technology and precise workmanship.

Advantages

F4 mobile driving system

The table has rubber adjustable feet and four casters for easy movement. To adjust the legs and position the table on the casters, use the foot lever.

Electric height adjustment

The table's height can be adjusted using a frame located in the lower part. This feature can be operated with the foot from any side of the table, providing convenience for the physiotherapist.

High-quality upholstery

The Theraspace tables feature high-quality upholstery that comes in a variety of colors. The upholstery is easy to clean and resistant to scratches, scrapes, and liquids such as blood, sweat, and urine.

Remote control

The table's height can be controlled using a wireless handheld remote control. To adjust the table's height, press the corresponding button on the remote control.

Position memory feature

The table has rubber adjustable feet and four casters for easy movement. To adjust the legs and position the table on the casters, use the foot lever.

Position locking system

The table is initially set to lock the movement of the actuators. If no movement is initiated after unlocking, the system will revert back to the locked state after ten seconds.

Lux upholstery

Choose the color most suitable for your office. The upholstery material is waterproof, antibacterial and easy to clean.

Temperature resistant

The two-layer upholstery of the **Theraspace Z3** massage table is resistant to high temperatures.

Liquid resistant

The top material of the **Theraspace Z3** massage table is resistant to liquids (water, blood, brine).

Comfortable surface

The **Theraspace Z3** massage table is covered with a soft, flexible and comfortable coating.

Scratch resistant

The two-layer upholstery of the **Theraspace Z3** features a scratch-resistant structure.

Accessories

Wedges

Bolsters

Paper holder

Remote controller

Stabilization Belt

Video presentation:

Product codes:

Reference: -

EAN13: -

UPC: -

Product features:

Type of table: Rehabilitation table

Number of sections: 3

Wheeling system: F4

Pivot: Yes

Length [cm]: 205

Width [cm]: 67

Height [cm]: 49-95

Weight [kg]: 110

Max. lifting capacity [kg]: 200

Max. working load [kg]: 275

Product attributes:

Tables - Color: OM006 - White, OM240 - Beige, OM315 - Blue, OM221 - Dark blue, OM314 - Grey, OM076 - Graphite, OM173 - Brown, OM225 - Black, OM469 - Light green, OM157 - Green, OM053 - Red, OM217 - Burgundy, OM095 - Mostaza, OM075 - Orange

Tables - Paper towel holder: No, Yes

Tables - Heated table surface: No, Yes

Tables - Remote Control: No, Yes

Power supply [V / Hz]: 230 / 60

Warranty: 36 months (Upholstery included)

Medical device: Yes